


# Locket

Suture Retention Device

A cost-effective solution for access site closure, ensuring efficient hemostasis while prioritizing patient comfort.



Allows for Same Day Hospital Discharge



Close multiple adjacent access sites with one device



Suitable for a wide range of catheter sizes up to 27F



Patient Comfort




Cost Effective




100% Immediate Hemostasis

# LockeT


Intuitive Suture Retention Device



**Step 1: Create a standard figure of eight and remove sheaths.**



**Step 2: Tighten sutures, cleanse site and place sutures within LockeT.**



**Step 3: Rotate the handle 90 degrees.**

## ORDERING INFORMATION

MODEL NUMBER	PRODUCT DESCRIPTION	QUANTITY
9301	Suture Retention Device	1 Box (10 devices/box)

LockeT is a Class 1 sterile product registered with the FDA under Device ID D496327

Refer to IFU for complete Indications for Use, contraindications, warnings, and precautions.



PN 82-9300-001 Rev03  
©2024 Catheter Precision, Inc. All Rights Reserved

 Catheter  
Precision

1670 Highway 160 West | Suite 205 | Fort Mill, SC 29708 USA  
T: (+1) 973-691-2000 | F: 973-691-7573  
[www.CatheterPrecision.com](http://www.CatheterPrecision.com)